

# Podiatry

---

["PBS" Diamond Burs](#)

---

[Accu-Mold™ - Moldable Silicone Compound](#)


---

[Medicaments](#)

---

[Nail Instruments](#)

---



## “PBS” Diamond Burs - Nail Reduction & Tissue Debridement

Professional Treatment System For Hyperkeratotic Tissue And Hypertrophic Nails. Through the exclusive bonding process, which permanently bonds diamond to the bur surface, Premier burs offer unsurpassed product life and performance! More efficient than steel or carbide instruments and with greater patient acceptance because of reduced heat generation.

All diamond burs fitted with  $\frac{3}{32}$ " shaft size.



**Barrel** • 2026005  
(To reduce thick, deformed nails.)



**Taper** • 2026030  
(For the treatment of incurvated nails or the nail groove.)



**Pear** • 2026010  
(To smooth the rough edges of nails.)



**Umbrella** • 2026035  
(For the reduction of hyperkeratoses, tylomas and heel callus.)




**Bud** • 2026015  
(The most popular pattern for nail reduction.)



**Small Olive** • 2026040  
(For working smaller areas.)



**Short** • 2026020  
(For better control in treating smaller digits.)



**Small Cylinder** • 2026045  
(For gross reduction of large nail areas.)



**Olive** • 2026025  
(Similar to the PEAR and BUD.)



**Diamond Bur System Kit** • 2026050  
A convenient kit of four nail and callus skin burs plus technique information. Includes one each:  
Bud, Short, Taper, Umbrella



**Clean-A-Diamond** • 2014001  
Dressing stone for diamond instruments.

Contact Premier Customer Service for bur cleaning and care guidelines.

## Accu-Mold™ - Moldable Silicone Compound



9045050 Accu-Mold - 200 Gram Kit (yields approximately 15 devices)

Legendary in its convenience, Accu-Mold™ remains a significant advance over other moldable compounds. Simple to mix, easy to mold, Accu-Mold™ forms to a firm yet pliable digital appliance with no mess or odor. Carefully following the detailed product instructions will produce a perfect digital appliance in just two to three minutes.

Accu-Mold™ has the following characteristics:

- Hypoallergenic, silicone based material.
- Once formed, retains original size and shape.
- Device can be cleansed with soap and water.
- Easily reshaped by grinding or trimming.

Accu-Mold™ is guaranteed for a shelf life of two years.

## Medicaments



Clinically proven and popular with Podiatrists, Lumicain® is a colorless hemostatic agent that provides fast and effective control of minor hemorrhage and granulation tissue formation.

9045010 Lumicain®



For the safe and effective treatment of verruca, Salinocaine contains a combination of salicylic acid and an anesthetic for patient comfort. Salinocaine is approved for topical use only.


9045015 Salinocaine

# Nail Instruments


Made In England



**Nail Nipper, SS, 5 1/2"**  
• 9085323 (Straight Jaw)



**Nail Nipper, SS, 5 1/2"**  
• 9085324 (Curved Jaw)



**Anvil Nail Split Forceps, SS, 5 1/2"**  
• 9085334

**Littauer Cutting Forceps, SS, 4"** • 9085393

**Littauer Cutting Forcep, SS, 5"** • 9085395

**Littauer Cutting Forceps, SS, 6"** • 9085397



# Nail Instruments

Made In Germany

**Collar-Crown 5"** • 9085346 (Curved/Sharp)

**Collar-Crown 5"** • 9085347 (Curved/Blunt)

Curved/Sharp Shown



**Collar-Crown 4"**

• 9085344 (Curved/Sharp)



**Collar-Crown 5"** • 9085342 (Straight/Sharp)



**Collar-Crown 4"**

• 9085345 (Curved/Blunt)



**Collar-Crown 4"** • 9085340 (Straight/Sharp)

**Collar-Crown 4"** • 9085341 (Straight/Blunt)

Straight/Sharp Shown

